

# Random Microseed Matrix Screening (rMMS): A new technique where the method is applied to membrane proteins in Lipidic Cubic Phase

Stefan Andrew Kolek<sup>a</sup>, Bastian Bräuning<sup>b</sup>, and Patrick Douglas Shaw Stewart<sup>a</sup>.


<sup>a</sup> Douglas Instrument Ltd, East Garston, Hungerford, Berkshire, RG17 7HD, United Kingdom

<sup>b</sup>Department of Chemistry, Technische Universität München, Lichtenbergstrasse 4, Garching, 85748, Germany.


## Introduction

- Inspired by the success of random Microseed Matrix Screening (rMMS), we have adapted rMMS to the crystallization of membrane proteins in LCP.
- LCP seed stock is made by scaling up LCP crystallization conditions without changing critical parameters.
- Seed crystals are grown directly in LCP, and (as with conventional rMMS) seeding is combined with an additive experiment.
- We used the method with the bacterial integral membrane protein *OmpF*: without microseeding, one new hit was found, but with LCP-rMMS eight new hits were found.
- We also demonstrate a method of generating seed gradients, which allows the number of crystals to be varied.


## Preparing the seed stock




## 1. Identify Initial Hit




- Screening experiments were dispensed to glass sandwich plates.
- We identified wells with high nucleation to make the seed stock.
- We used a published crystallization condition: 1.25 M K thiocyanate, 2.1 M Li nitrate and 0.1 M Na acetate pH 4.6. (Efremov & Sazanov, 2012)



## 2. Make the LCP Seed Stock



- Inject LCP into a syringe containing the crystallization cocktail.
- When the first sign of nucleation is observed, collect the LCP seed stock mixture and dispense straight away.
- It is important to collect the LCP seed stock before the crystals grow too much to ensure there is an crystallized protein in the LCP seed stock mixture.
- Using a 250 $\mu$ L or 500 $\mu$ L syringe allows more seed stock to be produced. It also improves reproducibility by better re-creating the crystallization environment and dimensions important to crystallization of a sandwich plate.




## 3. Dispense LCP rMMS Experiment

- Dispense the LCP seed stock to new screen.
- We found new hit conditions using the LCP-rMMS method.
- The main precipitant of the new hits were in all cases organic (PEG, ethanol or MPD), whereas the seed crystals were high-salt conditions.


Summary of new hits obtained in 120 conditions using LCP-rMMS.

Crystallization screen	Condition No.	Chemical composition	Without seeding		With seeding	
			Crystals?	Crystal size ( $\mu$ m)	Crystals?	Crystal size ( $\mu$ m)
JCSG+	17	40% MPD, 5% PEG 8K, 0.1 M sodium cacodylate pH 6.5	—	—	Yes	20–30
JCSG+	18	40% ethanol, 5% PEG 1K, 0.1 M phosphate-citrate pH 4.2	—	—	Yes	20–30
JCSG+	22	50% PEG 200, 0.2 M MgCl <sub>2</sub> , 0.1 M sodium cacodylate pH 6.5	—	—	Yes	5–10
JCSG+	30	40% PEG 300, 0.1 M phosphate-citrate pH 4.2	—	—	Yes	10–20
JCSG+	43	40% PEG 400, 0.2 M lithium sulfate, 0.1 M Tris-HCl pH 8.5	—	—	Yes	10–20
JCSG+	53	40% MPD, 0.1 M CAPS pH 10.5	—	—	Yes	20–30
JCSG+	64	20% Jeffamine M-600 pH 7.0, 0.1 M Na HEPES pH 7.5	Yes	10–20	—	—
JCSG+	66	10% MPD, 0.1 M Na bicine pH 9.0	—	—	Yes	10–20
Membrane 1	5	48% PEG 400, 0.2 M CaCl <sub>2</sub> , 0.1 M HEPES pH 7.5	—	—	Yes	10–20



## Preparation of LCP Gradient

- The ability to control the number of crystals per drop by diluting the seed stock is an important advantage of the rMMS method as applied to soluble proteins.
- An LCP seed stock gradient can be produced by mixing LCP seed stock and LCP without seeds together into a single syringe.
- The mixture was dispensed into a sandwich plate where all wells contained an aqueous solution that had previously given crystals with the LCP seed stock.
- A gradient effect was observed, however the length of the gradient was comparatively short, resulting in wells with clusters of crystals within clear LCP.



## Key Points

- We conclude that an LCP seed stock can easily be produced within gas-tight syringes by following a simple protocol.
- This protocol allows LCP crystallization conditions to be scaled up very easily and reliably because the physical dimensions that are critical for crystallization are maintained.
- We tested the method with the membrane protein *OmpF* from *E. coli*. The method worked well, increasing the number of hits from one to nine and identifying several new precipitants that supported crystallization.
- It is possible to dispense a LCP seed stock gradient, however we would have preferred to see individual crystals that were well separated from other crystals to encourage the growth of large crystals.

