

# Understanding “random” MMS microseeding

## – how new strategies can improve productivity

(A review of published work including our own)

Patrick Shaw Stewart, Stefan Kolek, Richard Briggs, Peter Baldock, Douglas Instruments Limited, near Oxford, UK:

### Introduction to “random” microseeding

#### Protein crystallization

- Step 1: screening with random solutions that have given crystals before x 96
- Step 2: optimization by making small changes
- Step 1.5: random microseeding

Modify your protein or make a new construct

Microseeding slide 1

#### Case study – Obmolova et al, Acta Cryst (2010) D66, 927 - 933

**Conventional methods**

Complexes: IL-13/C836 (mouse antibody) **No hits**

IL-13/H2L6 (humanized mAb) **One hit**

IL-13/M1295 (affinity-matured humanized mAb) **No hits**

**Random microseeding (RMS)**



Optimization: 2.0 Å res. monoclinic P2<sub>1</sub>

Optimization: Both 1.9 Å resolution orthorhombic P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub>

Optimization: 2.8 Å res. P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub>

Optimization: 2.8 Å res. P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub>

#### Microseeding can be used



### CASE STUDY

Obmolova, Galina, et al. "Protein crystallization with microseed matrix screening: application to human germline antibody Fabs." *Acta Cryst. F* 70.8 (2014): 1107-1115.

#### Cross-seeding gave better diffraction

Cross-seeding is very helpful

VH/VL	B3	A27	L6	O12
1-69	○	○	○	○
3-23	○	○	○	○
3-53	○	○	○	○
5-51	○	○	○	○

Initial hit: PEG and AS

Self-seeding into PEG and AS

2.8 Å


rMMS with cross-seeding into PEG and Li sulfate

2.3 Å!

Microseeding slide 11



#### Strategy: mix crystal seeds from homologous proteins

VH/VL	B3	A27	L6	O12
1-69	○	○	○	○
3-23	○	○	○	○
3-53	○	○	○	○
5-51	○	○	○	○




Microseeding slide 12

### Microseeding with LCP



The seed conditions: OmpF in LCP with 1.1 M KSCN, 2.1 M LiNO<sub>3</sub>, 0.1 NaOAc pH 4.6

Screen	Cond	Main precipitant	Count	Comment	Without seeding	With seeding
JCSG+	17	40% MPD	0		0	1 Small crystals
JCSG+	18	40% ethanol	0		0	1 Small crystals
JCSG+	22	50% PEG 200	0		0	1 Very small crystals
JCSG+	30	40% PEG 300	0		0	1 Very small crystals
JCSG+	43	40% PEG 400	0		0	1 Very small crystals
JCSG+	53	40% MPD	0		0	1 Small crystals
JCSG+	64	20% Jeffamine	1		0	1 Very small crystals
JCSG+	66	10% MPD	0		0	1 Very small crystals
Memb 1	5	48% PEG 400	0		0	1 Very small crystals
Total			1		0	8



Results: new crystallization conditions found.


Top: no seed stock added

Bottom: with seed stock

### Crystallizing complexes – cross-seed with crystals of one component


#### Cross-seeding

You don't have to match the unit cell, only one of the structural planes of the crystals



Microseeding slide 19


#### Focusing on “receptive” conditions



Microseeding slide 22

#### Crystallizing complexes

Radaev and Sun. Crystallization of protein-protein complexes. *J. Appl. Cryst.* (2002). 35, 674-676



PEG / (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> / other salts / organic solvents (including 2-propanol, MPD, ethanol)

Random uncomplexed samples, all protein-protein complexes included in this survey, immune complexes, antibody-antigen complexes, signal transduction complexes, receptor and ligand complexes, miscellaneous protein-protein complexes, enzyme related complexes, oligomeric protein complexes

Microseeding slide 21

#### Appearance of uncrushed crystals after incubation for one day

Protein	Crystals in Hit Sol.	Crystals in Isopropanol	Crystals in PEG 600	Crystals in Amm. sul.	Crystals in NaCl	Crystals in protein stock
Gluc Isom	OK	Cracked	Shattered	Cracked	Dissolved	Dissolved
Hemoglobin	OK	Cracked	OK	Dissolved	Dissolved	Dissolved
Thaumatin	OK	Cracked	OK	OK	OK	Grew
Thermolysin	OK	OK	Shattered	OK	Dissolved	Grew
Trypsin	OK	OK	Dissolved	OK	OK	Dissolved
Xylanase	OK	OK	Cracked	OK	OK	Dissolved


Microseeding slide 26

### Experimental design for finding the right dilution of the seed stock

- If solutions to be optimized are arranged in columns up to 8 hits can be optimized in one experiment/plate
- 10<sup>-3</sup> to 10<sup>-6</sup> dilutions were required for 3 model proteins

- “Additive Scatter” script can be used with 2-d grids
- Seed stocks vary orthogonally
- Up to 5 seed dilutions per plate
- Any size block can be filled


#### New “combinatorial” experimental design




#### Microseeding:

- A1: 100% seed stock
- A2: 25% seed stock
- A3: 6.3% seed stock
- A4: 1.6% seed stock
- A5: 0.4% seed stock
- A6: 0.1% seed stock
- A7: 0.02% seed stock
- A8: 0.006% seed stock
- A9: 0% seed stock (control)

#### New “combinatorial” experimental design



Microseeding slide 32



Microseeding slide 31